


Launching Simple Viewer

For this example guide all locations will be the "Test" folder as shown below in the Launching Simple Viewer section.

- Launch the **Simple Viewer Launcher** as shown below.



Choosing Viewer Type

Once the Simple Viewer Launcher opens you will have the choices below.

Opening the **Simple Viewer Original** will open the native Ez3D-i viewer software.

Note: Graphics card requirements are a NVidia GeForce GT 720 or higher to use the Simple Viewer Original. If your computer's graphics card does not meet these specifications, use the Simple Viewer Lite.

✕

Ez3D-i Simple Viewer

Simple Viewer Original requires a dedicated graphics card (NVIDIA GeForce GT 720 or above).
Simple Viewer Lite does not require a dedicated graphics card, but work done before exporting the image from Ez3D-i cannot be viewed with Simple Viewer Lite.

Simple Viewer Original

Simple Viewer Lite

ewoosoft

[Read Me](#)