

1. General Information

1.1. Equipment Overview

1.1.1. Indications for Use

Green Smart (Model: PHT-35LHS) is intended to produce panoramic, cephalometric or 3D digital x-ray images. It provides diagnostic details of the dento-maxillofacial, sinus and TMJ for adult and pediatric patients. The system also utilizes carpal images for orthodontic treatment. The device is to be operated by healthcare professionals.

1.1.2. Main Features

Green Smart is an advanced 4-in-1 digital X-ray imaging system that incorporates PANO, CEPH (Optional), CBCT and 3D MODEL Scan (Optional) imaging capabilities into a single system and its main features are as follows.

- Multi FOV support: Selectable FOV among 10x8.5, 10x7, (Optional) 5x5 (cm)
- Multi imaging solution for accurate diagnostics
- Conventional 2D (PANO, CEPH) image acquisition
- 3D scanning for Plaster Cast with FOV 10x8.5 (cm) (Optional)
- DICOM (Digital Imaging Communication in Medicine) format supported

1.2 Available Modes

Equipment	Available Modes
Green Smart	PANO, CEPH (Optional), CBCT, 3D MODEL Scan (Optional)

2. Functional Specifications

2.1. PANO Mode

2.1.1. Overview

PANO imaging software is classified into two levels as below.

Level	Examination Option	Optional Status	
		Domestic	Overseas
Normal	PANO Examination	Default	Default
	Special Examination	Default	Default
Magic PAN	Applies to entire PANO Examination programs	Default	Optional

2.1.2. Scan Time / Exposure Time

Examination Mode	Arch Type	ROI	High Resolution		Green (Optional)	
			Scan Time (s)	Exposure Time (s)	Scan Time (s)	Exposure Time (s)
PANO Examination	Narrow	Standard	14.1	13.5	7.0	7.0
		Right	14.1	6.8	7.0	3.5
		Front	14.1	11.3	7.0	5.8
		Left	14.1	6.8	7.0	3.5
	Normal	Standard	14.1	13.5	7.0	7.0
		Right	14.1	6.8	7.0	3.5
		Front	14.1	11.3	7.0	5.8
		Left	14.1	6.8	7.0	3.5
	Wide	Standard	14.1	13.5	7.0	7.0
		Right	14.1	6.8	7.0	3.5
		Front	14.1	11.3	7.0	5.8
		Left	14.1	6.8	7.0	3.5
	Child	Standard	12.0	11.5	6.8	6.7
		Right	12.0	5.7	6.8	3.3
		Front	12.0	9.2	6.8	5.2
		Left	12.0	5.7	6.8	3.3

Examination Mode	Arch Type	ROI	High Resolution		Green (Optional)	
			Scan Time (s)	Exposure Time (s)	Scan Time (s)	Exposure Time (s)
PANO Examination	Orthogonal	Standard	14.1	13.5	7.0	7.0
		Right	14.1	6.7	7.0	3.5
		Front	14.1	11.1	7.0	5.7
		Left	14.1	6.7	7.0	3.5
		Bitewing	14.1	9.2	7.0	5.0
		Bitewing Incisor (Optional)	14.1	2.8	7.0	1.4
		Bitewing Right	14.1	5.0	7.0	2.8
		Bitewing Left	14.1	5.0	7.0	2.8
SPECIAL Examination	-	TMJ LAT Open	14.1	6.7	14.1	6.7
		TMJ LAT Close				
		TMJ PA Open (Optional)	10.0	6.1	10.0	6.1
		TMJ PA Close (Optional)				
		Sinus LAT (Optional)	4.0	3.7	4	3.7
		Sinus PA	8.8	7.7	8.8	7.7


- *Scan Time: The actual time that the equipment shoots the patient except for the initial acceleration and late deceleration stages.*
- *Exposure Time: The actual time that the patient is exposed to the X-ray emission*

2.1.3. Exposure Condition

Examination Mode	Image Quality	Gender / Age group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
PANO Examination	High Resolution	Man	Hard	75	12
			Normal	74	12
			Soft	73	12
		Woman	Hard	74	12
			Normal	73	12
			Soft	72	12
		Child	Hard	68	10
			Normal	67	10
			Soft	66	10
	Green (Optional)	Man	Hard	75	14
			Normal	74	14
			Soft	73	14
		Woman	Hard	74	14
			Normal	73	14
			Soft	72	14
Child		Hard	68	12	
		Normal	67	12	
		Soft	66	12	
SPECIAL Examination	N/A	Man	Hard	75	14
			Normal	74	14
			Soft	73	14
		Woman	Hard	74	14
			Normal	73	14
			Soft	72	14
		Child	Hard	68	12
			Normal	67	12
			Soft	66	12

- Recommended exposure condition can be different from the values applied to the equipment.

2.1.4. Sample Images

Item	Images
PANO Examination	 <p data-bbox="790 846 1125 880"><High Resolution/ Standard></p>

2.2. CEPH Mode

2.2.1. Scan Time / Exposure Time

Examination Mode	High Resolution		Green	
	Scan Time (s)	Exposure Time (s)	Scan Time (s)	Exposure Time (s)
Lateral	3.9	3.9	1.9	1.9
Full Lateral	5.4	5.4	3.9	3.9
PA	4.9	4.9	2.4	2.4
SMV	4.9	4.9	2.4	2.4
Waters' view	4.9	4.9	2.4	2.4
Carpus	4.9	4.9	2.4	2.4

- *Scan Time: The actual time that the equipment shoots the patient except for the initial acceleration and late deceleration stages.*
- *Exposure Time: The actual time that the patient is exposed to the X-ray emission.*



2.2.2. Exposure Condition

Examination Program	Image Option	Gender / Age group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
Lateral	High Resolution)	Man	Hard	92	15.0
			Normal	90	15.0
			Soft	88	15.0
		Woman	Hard	90	15.0
			Normal	88	15.0
			Soft	86	15.0
		Child	Hard	88	15.0
			Normal	86	15.0
			Soft	84	15.0
	Green	Man	Hard	92	16.0
			Normal	90	16.0
			Soft	88	16.0
		Woman	Hard	90	16.0
			Normal	88	16.0
			Soft	86	16.0
Child	Hard	88	16.0		
	Normal	86	16.0		
	Soft	84	16.0		

Examination Program	Image Option	Gender / Age group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
Full Lateral	High Resolution / Green	Man	Hard	92	14.0
			Normal	90	14.0
			Soft	88	14.0
		Woman	Hard	90	14.0
			Normal	88	14.0
			Soft	86	14.0
		Child	Hard	88	14.0
			Normal	86	14.0
			Soft	84	14.0
PA SMV Waters' view	High Resolution	Man	Hard	92	14.0
			Normal	90	14.0
			Soft	88	14.0
		Woman	Hard	90	14.0
			Normal	88	14.0
			Soft	86	14.0
		Child	Hard	88	14.0
			Normal	86	14.0
			Soft	84	14.0
	Green	Man	Hard	92	15.0
			Normal	90	15.0
			Soft	88	15.0
		Woman	Hard	90	15.0
			Normal	88	15.0
			Soft	86	15.0
Child		Hard	88	15.0	
		Normal	86	15.0	
		Soft	84	15.0	
Carpus	High Resolution / Green	Man	Hard	90	6.0
			Normal	88	6.0
			Soft	86	6.0
		Woman	Hard	88	6.0
			Normal	86	6.0
			Soft	84	6.0
		Child	Hard	86	6.0
			Normal	84	6.0
			Soft	82	6.0

● Recommended exposure condition can be different from the values applied to the equipment.

2.2.3. Sample Images

Item	Images
Lateral	 <p data-bbox="911 1061 1002 1088"><Green></p>
PA	 <p data-bbox="911 1792 1002 1818"><Green></p>

2.3. CBCT Mode

2.3.1. Exposure Area

FOV (cm)	Vertical Position	Horizontal Position		
		Right	Center	Left
10x8.5	Occlusion	X	O	X
10x7	Occlusion	X	O	X
5x5	Maxilla / Mandible	Right Molar / Right / Incisor / Left / Left Molar		

- FOV 10x8.5 and 10x7 capture full arch area. (10x8.5 : Adult mode, 10x7 : Child mode)

2.3.2. Scan Time / Exposure Time

FOV (cm)	Scan Time (s) (High Resolution / Green)	Exposure Time (s) (High Resolution / Green)
10x8.5	18.0	16.4
10x7	18.0	16.4
5x5	18.0	11.0

- Scan Time: The actual time that the equipment shoots the patient except for the initial acceleration and late deceleration stages.
- Exposure Time: The actual time that the patient is exposed to the X-ray emission.

2.3.3. Exposure Condition

FOV (cm)	Image Quality	Gender / Age Group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
10x8.5	High Resolution	Man	Hard	95	8.7
			Normal	94	8.7
			Soft	93	8.7
		Woman	Hard	95	8.4
			Normal	94	8.4
			Soft	93	8.4
		Child	Hard	95	8.1
			Normal	94	8.1
			Soft	93	8.1
	Green	Man	Hard	80	5.0
			Normal	79	5.0
			Soft	78	5.0
Woman		Hard	80	4.7	
		Normal	79	4.7	
		Soft	78	4.7	

FOV (cm)	Image Quality	Gender / Age Group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
10x7	High Resolution	Child	Hard	80	4.4
			Normal	79	4.4
			Soft	78	4.4
		Man	Hard	95	8.7
			Normal	94	8.7
			Soft	93	8.7
		Woman	Hard	95	8.4
			Normal	94	8.4
			Soft	93	8.4
	Child	Hard	95	8.1	
		Normal	94	8.1	
		Soft	93	8.1	
	Green	Man	Hard	80	5.0
			Normal	79	5.0
			Soft	78	5.0
Woman		Hard	80	4.7	
		Normal	79	4.7	
		Soft	78	4.7	
Child		Hard	80	4.4	
		Normal	79	4.4	
		Soft	78	4.4	
5X5	High Resolution	Man	Hard	95	8.7
			Normal	94	8.7
			Soft	93	8.7
		Woman	Hard	95	8.4
			Normal	94	8.4
			Soft	93	8.4
	Child	Hard	95	8.1	
		Normal	94	8.1	
		Soft	93	8.1	
Green	Man	Hard	90	6.2	
		Normal	89	6.2	
		Soft	88	6.2	
	Woman	Hard	90	5.9	
Normal		89	5.9		



FOV (cm)	Image Quality	Gender / Age Group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
			Soft	88	5.9
		Child	Hard	90	5.6
			Normal	89	5.6
			Soft	88	5.6

2.3.4. Reconstruction Time / File Size

FOV(cm)	Voxel Size (mm)	Reconstruction Time (s)	File Size (MB)
10x8.5	0.2	103	260
	0.3	107	77
10x8.5 (with Auto Pano*)	0.2	106	267
	0.3	109	84
10x7	0.2	85	202
	0.3	85	60
5x5	0.08	123	466
	0.12	119	139


- The above data is obtained from a computer system which is based on HP Workstation Z440 and Geforce GTX1060 6GB
- Image reconstruction time varies depending on computer specifications and/or working conditions.

2.3.5. Sample Images

FOV (cm)	Images
10X8.5	 <p data-bbox="799 976 1133 1010"><High Resolution, Voxel 0.2></p>
10X7	 <p data-bbox="799 1621 1133 1655"><High Resolution, Voxel 0.2></p>

FOV (cm)	Images
----------	--------


5X5



< Mandible – Right Molar : High Resolution, Voxel 0.08 >



< Mandible – Right : High Resolution, Voxel 0.08 >

FOV (cm)	Images
5X5	 <p data-bbox="671 891 1259 929">< Mandible – Incisor : High Resolution, Voxel 0.08 ></p>

대학원

2.4. 3D MODEL Scan Mode

2.4.1. Exposure Area

FOV (cm)	MODEL Type	Horizontal Position		
		Right	Center	Left
10x8.5	Upper (Maxilla)	X	O	X
	Lower (Mandible)	X	O	X

2.4.2. Scan Time / Exposure Time

FOV (cm)	Scan Time (s)	Exposure Time (s)
10x8.5	18.0	16.4

- *Scan Time: The actual time that the equipment shoots the patient except for the initial acceleration and late deceleration stages.*
- *Exposure Time: The actual time that the patient is exposed to the X-ray emission.*

2.4.3. Exposure Condition

FOV (cm)	Gender / Age Group	X-ray Intensity	Tube Voltage (kVp)	Tube Current (mA)
10x8.5	Man / Woman / Child	Hard / Normal / Soft	95	8.7

- *Recommended exposure condition can be different from the values applied to the equipment*

2.4.4. Reconstruction Time / File Size

FOV(cm)	Voxel Size (mm)	Reconstruction Time (s)	File Size (MB)
10x8.5	0.12	179	1203

- *HP Workstation Z440 and NVIDIA Geforce GTX1060 6GB*
- *Image reconstruction time varies depending on computer specifications and/or working conditions.*



NOTE

The 3D MODEL Scan modality is not available for EasyDent / Ez3D Plus users.

3. PC Specification (Recommended)



Item	Specifications (HP)
CPU	Intel Xeon E5-1607v3 4C 3.1GHz 1866 4C or Higher
RAM	2X8GB DDR4-2400 Registered RAM
Hard disk drive	1TB SATA 1 st HDD
Graphic board	GALAXY Geforce GTX1060 6GB
Ethernet interface	10/100/1000 Mbps, RJ-45, 2 Port
Serial Port (RS232)	HP Serial Port Adapter Kit
Power supply	≥ 700 Watts (90% Efficiency)
Slots	1 PCI Express Gen3 x16 slot 1 PCI Express Gen3 x8 slot 1 PCI Express Gen2 x4 Slot 1 PCI Express Gen2 x1 Slot 1 PCI Slot
CD/DVD drive	DVD-ROM, DVD+/-RW, Blu-Ray
Monitor Resolution	1280 x 1024 screen resolution
Operating system	Windows 10 Professional 64-Bit OS
Bios	2.14
Recommended system	Z440

- *Recommended PC specification can be changed without prior notice.*

4. Mechanical Specification

4.1. Equipment Dimensions

[Unit: mm]



Item		Description
Weight	Without CEPH unit	143 kg (315.3 lbs. – without base)
		196 kg (432.1 lbs. – with base)
	With CEPH unit	168 kg (370.4 lbs. - without base)
		221 kg (487.2 lbs. - with base)
Total Height	Without base	Max. 2304.5 mm
	With base	Max. 2336mm
Dimensions during operation (Length x Width x Height)	Without CEPH unit	1389.3 mm (L) x 1113 mm (W) x 2304.5 mm (H) (without base)
		1389.3 mm (L) x 1113 mm (W) x 2336 mm (H) (with base)
	With CEPH unit	1390.3 mm (L) x 1882.3 mm (W) x 2304.5 mm (H) (without base)
		1390.3 mm (L) x 1882.3 mm (W) x 2336 mm (H) (with base)

Item	Description
Rotating Unit Vertical Movement	Max. 800 mm
Installation Type	Base Stand / Wall Mount (Default: Wall Mount type)
Packing Box Organization	Main Box, CEPH Box (Optional), Base Box (Optional)

4.2. Image Magnification

Mode	FDD (mm)	FOD (mm)	ODD (mm)	Magnification
PANO	600	477.7	122.3	1 : 1.25
CEPH	1745	1524	221	1 : 1.14
CBCT	600	428.6	171.4	1 : 1.4

* FDD : Focal Spot to Detector Distance

* FOD : Focal Spot to Object Distance

* ODD : Object to Detector Distance ($ODD = FDD - FOD$)

* Magnification = FDD / FOD

5. Technical Specification

5.1. X-ray Generator

Item		Description	
Generator Specification	Model	DG-07E22T2	
	Rated output power	1.6 kW	
	Inverter model name	INV-22	
	Type	Inverter	
	Normal/ Pulse	kVp	60 kV ~ 99 kV (1 kV increment)
		mA	4 mA ~ 16 mA (0.1 mA increment for CBCT, 1 mA increment for PANO and CEPH)
	Cooling	Automatically controlled / Protect $\geq 60^{\circ}\text{C}$ (Cooling fan: Optional)	
	Total filtration	Min. 2.5 mm Al	
	Added filter	1.5 mm Al (Fixed) / PANO and CEPH mode 1.5 mm Al (Fixed) + 3.0 mm Al (Automatically added) / CBCT mode	
	X-ray Tube	Manufacturer	Toshiba
Model		D-052SB (Stationary Anode type)	
Focal spot size		0.5X0.5mm(IEC 60336)	
Target Angle		5 degree	
Inherent Filtration		At least 0.8 mm Al equivalent at 50 kV	
X- ray Coverage		95 mm x 380 mm at SID 550 mm	
Anode Heat Content		35 kJ	
Duty Cycle		1:60 or more (Exposure time : Interval time)	

● Test Condition

Mode	Tube Voltage (kVp)	Tube Current (mA)	Exposure Time (s)
PANO	60 ~ 90	4 ~ 14	13.5
	60 ~ 90	4 ~ 14	11.5
	60 ~ 90	4 ~ 14	11.3
	60 ~ 90	4 ~ 14	11.1
	60 ~ 90	4 ~ 14	9.2
	60 ~ 90	4 ~ 14	7.7
	60 ~ 90	4 ~ 14	7.0
	60 ~ 90	4 ~ 14	6.8
	60 ~ 90	4 ~ 14	6.7
	60 ~ 90	4 ~ 14	6.1
	60 ~ 90	4 ~ 14	5.8
	60 ~ 90	4 ~ 14	5.7
	60 ~ 90	4 ~ 14	5.2
	60 ~ 90	4 ~ 14	5.0
	60 ~ 90	4 ~ 14	3.7
	60 ~ 90	4 ~ 14	3.5
	60 ~ 90	4 ~ 14	3.3
	CEPH	60 ~ 99	4 ~ 14
60 ~ 99		4 ~ 14	1.4
60 ~ 99		4 ~ 16	1.9
60 ~ 99		4 ~ 15	2.4
60 ~ 99		4 ~ 15	3.9
CBCT	60 ~ 99	4 ~ 14	4.9
	60 ~ 99	4 ~ 14	5.4
	60 ~ 99	4 ~ 12	16.4
	60 ~ 99	4 ~ 12	11.0

5.2. Detector Specifications

Item	Description	
	PANO & CBCT	CEPH
Model	Xmaru1404CF-PLUS	Xmaru2602CF
Detector Type	CMOS photodiode array	
Pixel size	198 μm @ 4X4 binning (49.5 μm @ no binning)	200 μm @ 2x2 binning (100 μm @ no binning)
Active area	CBCT – 135.8 mm X 36.4 mm PANO – 135.8 mm X 5.9 mm	259mm x 15.6mm
Frame Rate	~308 (Maximum FPS by Full Size Mode(686x184 @ 4x4 binning))	~330 @ 2x2 binning
Analogue-Digital Conversion	14 bits	
Operating condition	10 ~ 35 °C (Temperature) / 10 ~ 75 % (Humidity)	
Storage condition	-10 ~ 60 °C (Temperature) / 10 ~ 75 % (Humidity)	
Sensor size	230(W) x 160(L) x 26(H) mm	279(W) x 110(L) x 20(H) mm
Sensor weight	1.5Kg	1.3Kg
Converter	CsI:Ti	
Energy Range	50 – 120 kVp	50 – 120 kVp
Readout	Charge amplifier array	
Video Output	Optic	
MTF	> 45% @ 1 lp/mm > 8% @ 2.5lp/mm	>35% @1 lp/mm >5% @ 2.5 lp/mm
DQE	> 60% @ ~0lp/mm	> 60% @ ~0lp/mm
Dynamic Range	> 80dB	\geq 70dB

6. Electrical Specifications

Item	Descriptions
Power supply voltage	100 - 240 V ~
Frequency	50/60 Hz
Power rating	2.0 kV
Accuracy	kVp \pm 10 %, Ma \pm 20 %, sec \pm (5 % + 50 ms)

- The input line voltage depends on the local electrical distribution system.
- Allowable input voltage fluctuation requirement : \pm 10%

7. Environmental Specification














Item	Descriptions	
During Operation	Temperature	10 ~ 35 °C
	Relative humidity	30 ~ 75 %
	Atmospheric pressure	860 ~ 1060 hPa
During Transport and Storage	Temperature	-10 ~ 60 °C
	Relative humidity	10 ~ 75 %
	Atmospheric pressure	860 ~ 1060 hPa

8. Software Information

Item	Model name	Interworking	Version for domestic market	Version for overseas market
Console Program	VCaptureSW	o	1.0.0.1	1.0.0.1
LCD Program	LCDSW	o	1.0.0.1	1.0.0.1
2D Viewer / Patient information management Program	EzDent-i	o	2.2	2.2
	EasyDent	x	4.5.2.7	4.1.5.9
3D Viewer	Ez3D-i	o	4.3	4.3
	Ez3D Plus	x	1.2.6.23	1.2.6.23

- SDK method has been applied for 3rd party programs.

9. Enclosed Components



Components	Name and Function	Materials
	Normal Bite : For PANO and CBCT normal patients	PC (Polycarbonate)
	Special Bite A : For PANO TMJ and Sinus modes & For CBCT TMJ patients	PC (Polycarbonate)
	Special Bite B : For PANO edentulous patients	PC (Polycarbonate)
	Normal Chinrest : For Normal Bite	ABS (Acrylonitrile butadiene styrene) copolymer
	Special Chinrest : For Special Bite A & Special Bite B	ABS (Acrylonitrile butadiene styrene) copolymer
	Temple Supports (1 set)	PC (Polycarbonate)
	Ear Rods (1 set)	Silicone
	Nasal Positioner Cover : For CEPH	Silicone
	Carpus Plate	PC (Polycarbonate)
	Sanitary Vinyl Covers (disposable) for the Bite 	LDPE (Low-density polyethylene)
	Protractor (1 set) : For positioning the patient's body in CEPH mode.	PC (Polycarbonate)
	Model Scan Jig	ABS (Acrylonitrile butadiene styrene) copolymer

10. Standards and Regulations

10.1. Standards

Green Smart is designed and developed to comply with the following international standards and regulations.

- IEC 60601-1, IEC/EN 60601-1-2, IEC 60601-1-3, IEC 60601-1-6, IEC 60601-2-63
- CAN/CSA-C22.2 No. 60601-1:14, CAN/CSA-C22.2 No. 60601-1-3:09, CAN/CSA-C22.2 No. 60601-1-6:11, CAN/CSA-C22.2 No. 60601-2-63:15, CAN/CSA-IEC 62366:15
- ANSI/AAMI ES60601-1:2005 / (R)2012, AND A1:2012, A2:2010 / (R)2012 (Consolidated text - edition 3.1)
- 21 CFR 1020.30, 31, 33
- NEMA Standard publication PS 3.1-3.18, 2008

	This is Class IIb equipment and obtained CE marking in April, 2007 for regulations compliance in accordance with the revised European Union's MDD (Medical Devices Directive) 93/42 EEC.
	This equipment received the CSA certification mark in accordance with CAN/CSA C22.2 No.601.1 regulations.

10.2. Classifications (IEC 60601-1 6.1)

- Degree of protection against water ingress: Ordinary Equipment: IPX0
- Degree of protection against electric shock: Class 1 equipment, Type B Applied Parts: Temple Supports, Chinrests and Bites.



11. Additional Information

For further information for Green Smart, please contact us at:

VATECH Co., Ltd.

Tel: +82-1588-9510

Website: www.vatech.co.kr

Headquarters: 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

Factory: 13, Samsung 1-ro 2-gil, Hwaseong-si, Gyeonggi-do, 18449, Korea

- ※ Due to a constant technological improvement, the contents of this Product Data may not contain the most updated information, subjecting to change without prior notice to the persons concerned.

VATECH RESERVES THE RIGHT TO CHANGE THE TERMS AND CONDITIONS OF THE DISTRIBUTION AND THESE TERMS AND CONDITIONS ARE SUBJECT TO CHANGE AT ANY TIME WITHOUT NOTICE. INFORMATION ON THIS PAGE IS THE CONFIDENTIAL PROPERTY OF VATECH GLOBAL. ANY UNAUTHORIZED DISTRIBUTION OR COPYING OF THIS PAGE IS STRICTLY PROHIBITED.